

HMP- SA400

STAND AID

HMP-SA400



The HMP-SA400 stand aid is a transport assist unit which keeps the resident active and engaged in the process. Transport functions of all types are quick and require minimal caregiver assistance. Users grasp the middle bar and pull themselves up into position. A padded split seat swings out to allow loading or unloading. The seat then swings back to form a secure and comfortable seat for toileting. An excellent alternative to wheelchairs for toileting. Rugged construction for years of service in acute care, nursing homes and private residences.



HMP-SB442 Optional Safety Strap

Made in Taiwan



Use with chair , commode

FEATURES:

- * Weight capacity: 400 lbs safety working load.
- * Dual knee pads for secure stand support and transports .
- * Sturdy yet easily fits through narrow doorways
- * Affordable alternative to electric stand aids.

Stand Aid Lift Specifications				
Maximum width	Maximum height	Overall length	Unit weight	Front Base
25"	42.5"	33.5"	61 Lbs.	5"

