

HMP- SL400/600

STAND ASSIST LIFT

HMP-SL400 Standard



FEATURES:

- * Weight capacity: 400/600 lbs safe working load.
- * Removable quick-release foot platform.
- * Superior comfort with deluxe padded slings.
- * Ergonomic steering handle for all caregivers
- * Foot pedal base opening








HMP-SL600 Bariatric



Made in Taiwan



SPECIFICATION:

WEIGHT CAPACITY	MAX. HEIGHT	MIN. HEIGHT	BASE CLOSED EXTERNAL WIDTH	BASE CLOSED INTERNAL WIDTH	BASE OPEN EXTERNAL WIDTH
 400/600 Lbs 400/600 Lbs (180/270 KGS)	 68" (172 CM)	 33" (84 CM)	 25.5" (65 CM)	 20.75" (53 CM)	 37" (94 CM)
BASE OPEN INTERNAL WIDTH	BASE HEIGHT CLEARANCE	MAST HEIGHT	OVERALL LENGTH	ACTUATOR	UNIT NET WEIGHT
 32.5" (83 CM)	 4.5" (11 CM)	 40" (102 CM)	 40" (102 CM)	 Linak # LA31/LA34	 117/122 Lbs (53/55 Kgs)

