

# Overbed Table



## HMP-S201

### FEATURES:

- ✘ W/ powder coating
- ✘ Height adj. from 710mm to 1175mm
- ✘ Freely position tilt wood panel
- ✘ 2 castor w/ brake
- ✘ 10kgs weight capacity on wood panel testing in 1175mm
- ✘ Wood panel dimension: 765x385x15mm



## HMP-S201W

### FEATURES:

- ✘ W/ wood panel with one piece injection PVC white color surface
- ✘ W/ powder coating
- ✘ Height adj. from 710mm to 1175mm
- ✘ Freely position tilt wood panel
- ✘ 2 castor w/ brake
- ✘ 10kgs weight capacity on wood panel testing in 1175mm
- ✘ Wood panel dimension: 600x400x18.50mm

### PACKAGING: HMP-S201/S201W

- ✘ 1pcs/ctn; Carton dimension: 645x70x650mm
- ✘ Container load: 864 ctns/1750ctns( 20foot/40foot)
- ✘ Weight: 8.5 kgs



## HMP-S202

### FEATURES:

- ✘ Wooden flame-resistant panel
- ✘ Chrome plated interior tube
- ✘ Powder coating exterior tube
- ✘ Portable castors with brakes
- ✘ Height adjustment: 715~1145 mm
- ✘ Weight capacity: 10 kgs



## HMP-S202W

### FEATURES:

- ✘ W/ wood panel with one piece injection PVC white color surface
- ✘ Chrome plated interior tube
- ✘ Powder coating exterior tube
- ✘ Portable castors with brakes
- ✘ Height adjustment: 718~1150 mm
- ✘ Weight capacity: 10 kgs



## HMP-S205

### FEATURES:

- ✘ Wooden flame-resistant panel
- ✘ Chrome plated interior tube
- ✘ Powder coating exterior tube
- ✘ Portable castors with brakes
- ✘ Height adjustment: 718~1150 mm
- ✘ Weight capacity: 10 kgs

